

Cryoport Express® Cryogenic HV3 Shipping System

The **Cryoport Express® Cryogenic HV3 Shipping System** has been reengineered to optimize storage capacity for air and ground transportation, operational efficiency, and mobility.

This liquid nitrogen dry vapor shipping system is available for general purpose shipments as well as larger cryogenic payloads and mitigates risk throughout the transportation process by leveraging our:

- Smartpak™ Condition Monitoring System
- Cryoport® Logistics Management System
- Live View™ Application
- 24/7/365 Customer Support

Recommended commodities: cell & gene therapies, vaccines, and other biologic materials.



Specifications

Temperature Range

-150°C or colder

Static Hold Time

Up to 13 days

Dynamic Hold Time

Up to 10 days

Exterior

22" L x 22" W x 28" H

Interior Chamber

8.5" diameter
12" height

Weight

Dry - 69 to 73 lb / 31 to 33 kg
Charged - 88 to 92 lb / 40 to 42 kg

Capacity - Canes

Each cane holds 5 vials

- 6 Safepak® bags - 15 canes per bag (75 cryovials per bag, 450 cryovials total)
- 100 canes (2.0 ml per vial)
- 4 metal canisters (500 cryovials total)
- Max packout can hold up to 500 cryovials

Capacity - Boxes

- 5 Standard cryovial boxes (5" x 5" x 2")
- 1 Safepak® XL with optional Cryostrap®
- 1 metal cryovial box rack (5 boxes per rack)

Capacity - Cassettes

- 9 cassettes (up to 750 ml); additional quantities depending on cassette size
- 1 Safepak® XL with optional Cryostrap®
- 1 metal cassette rack

Benefits

- With its retractable handle and robust wheels, the shipping system promotes safer maneuverability and may only require one operator as opposed to two.
- Mating top ridges and the formed base make it possible to properly stack and nest a shipping systems at storage facilities and sites of care.
- Modular 360° foam coverage within the enclosure provides biologic materials with additional protection, enhanced impact resistance, and reduced dewar movement.
- The contoured base design enhances the contact area between the enclosure foam and the dewar, thereby reducing impact pressure while improving the enclosure's structural strength.
- Fits and complies with narrow-bodied aircraft to minimize costs, optimize the required hold time, avoid shipment delays, reduce potential fees related to palletization, and prevent rejection due to enclosure design.

Contact Info

cryoport.com
info@cryoport.com
+1 949.232.1900



Stacking & Nesting Abilities

The interlocking top ridges and formed base allow for seamless stacking and nesting of two shipping systems, to maximize storage and available space during transportation. We recommend a maximum of two layers of stacking on a 48" x 48" pallet, which is equivalent to eight total shipping systems.

Additionally, the offset height ridges deter misorientation during both air and ground transport.



Key Features

Top Label Section: Airline airway bills; Client labels; 12.5" x 12.5" area

Integrated Handle: Front-facing for easy access mobility

Label Panel: Dedicated space for classification labels and markings

Double-action Steel Latches: Promotes security and positive engagement

Integrated Wheels: Robust and inset; easy mobility

Large Document Storage: Nests between 2 ribs and sits on the foam. Fits 10" x 12" plastic travel pouch for shipping documents and non-temperature sensitive accessories

Shipper Identification Label: All client labels featured on one side

Document Storage Dimensions: 13.5" D x 8.75" L x 5.75" W

