


CRYOPORT EXPRESS® HIGH VOLUME SLIDERITE® PALLET SYSTEM



General Use Instructions


NEED HELP? CONTACT CRYOPORT SYSTEMS' CLIENT CARE AT (949) 470-2305 OR EMAIL CS@CRYOPORT.COM

Product Overview



The Sliderite® is a mobile pallet that was designed to guarantee that the Cryoport Express® HV2 Shipping System remains in the upright position and is easily transported.

Modes



The Sliderite® pallet has two primary modes that are controlled by the pedals on either end of the pallet.

Pushing down on the **RED** pedal activates the stationary mode, which prevents the pallet from rolling.

Pushing on the **GREEN** pedal releases the pallet into motion mode, which allows it to roll around.