

# Shipment Order Form (GP)

*Cryoport Systems does NOT ship Category A commodities at this time.*

## Bill To:

Company

Address

City

State

Postal Code

Country

Account No.

Contact

Phone

Email

PO Number

Quote Number

## Ship From:

Company

Address

City

State

Postal Code

Country

Tax ID

Contact

Phone

Email

Expected Arrival Date and Time

## Pickup Information:

**Section A:** Schedule pickup. If not required, please go to Section B.

Pickup Date

Pickup Time

to

**Section B:** No Pickup required.

Does this location have a daily pickup?

Will the package be dropped off?

**Ship To:**

Company

Address

City

State

Postal Code

Country

Account No.

Contact

Phone

Email

Preferred Courier Type:

Standard Courier -  
FedEx, UPS, DHL (if possible)

Specialty Courier -  
If selected, please specify:

**International Shipments:**

Will use own broker:      Yes                  No

*If you selected Yes, please fill out our Broker Designation Form, Form-208*

**Shipping System & Quantity:**

For additional information and to download our specification sheets, please visit: [www.cryoport.com/brochures](http://www.cryoport.com/brochures)



**Cryoport Express® Cryogenic ST1 Shipping System**  
CXST1SPGP  
-150°C or colder



**Cryoport Express® Cryogenic CB1 Shipping System**  
CXCB1SPGP  
-150°C or colder  
\*The CB1 includes a top accessory tray for non-temperature sensitive items.



**Cryoport Express® Cryogenic HV3 Shipping System**  
CXHV3SPGP  
-150°C or colder



**Cryoport Express® C3™ Refrigerated Shipping System**  
C3CH2124SPGP  
2°C to 8°C



**Cryoport Elite® 28L Ultra Cold Shipping System**  
ELUC128SPGP  
-60°C to -80°C



**Cryoport Elite® 56L Ultra Cold Shipping System**  
ELUC156SPGP  
-60°C to -80°C

**Commodity Information:**

Description

Primary Container

If there are other secondary containers, please specify:

Container Volume

Quantity

Unit Value

Declared Value

Currency Code (On Declared Value)

**International Shipments:**

Country of Origin

Harmonized Code (Schedule B No.)

Reason for Shipment

If other, please specify:

**Commodity Information:**

1. **Exempt Human Specimen:** Specimen packages marked as "Exempt human specimen" according to the ICAO Technical Instructions are not regulated under the Hazardous Materials Regulations (HMR; 49 CFR Parts 171-180). In the United States, the mark "Exempt Human Specimen" is an indication that there is no infectious substance in the package. Packages bearing these marks may be accepted by an air carrier that has made a business decision not to accept hazardous materials.
2. **Exempt Animal Specimen:** Specimen packages marked as "Exempt animal specimen" according to the ICAO Technical Instructions are not regulated under the Hazardous Materials Regulations (HMR; 49 CFR Parts 171-180). In the United States, the mark "Exempt Animal Specimen" is an indication that there is no infectious substance in the package. Packages bearing these marks may be accepted by an air carrier that has made a business decision not to accept hazardous materials.
3. **Not Subject to Dangerous Good Regulations:** Substances which do not contain infectious substances or substances which are unlikely to cause disease in humans or animals and are not of human or animal origin.
4. **Genetically Modified Organism (GMO) (UN3245):** An organism whose genetic characteristics have been altered using the techniques of genetic engineering.
5. **Biological Substance, Category B (UN3373):** An infectious substance not in a form generally capable of causing permanent disability or life-threatening or fatal disease in otherwise healthy humans or animals when exposure to it occurs. This includes Category B infectious substances transported for diagnostic or investigational purposes.

***NOTE: Cryoport Systems does NOT ship Category A commodities at this time.***

## Cryoport Accessories

For additional information and to download our specification sheets, please visit: [www.cryoport.com/brochures](http://www.cryoport.com/brochures)

**Condition Monitoring** - Please check box to order

**Smartpak II®**  
ACC-9123

**Reporting** - Please check box to order

**Live View®**

**Live View® Lite - All Monitoring Data**  
SRV-10621

**Live View® Lite - Commodity Monitoring Data Only**  
SRV-10622

**Specimen Management** - Please insert quantity to order



**Cryovial Canes**

ACC-9033

Holds up to 5, 2.0 ml cryovials or equivalently sized alternative



**Cassette Rack**

ACC-9122

Holds 9 blood cassettes. Cassette sizes up to 11" H x 6" L x .5" D



**Cryovial Box Rack**

ACC-8905

Holds 5 cryovial boxes. Cryovial box sizes up to 5" W x 5" L x 2" H

**Packaging** - Please insert quantity to order



**Safepak®**

ACC-9139

Holds up to 15 cryovial canes. Holds 2, 65 mm goblets. Holds 2 Cryoport cryovial boxes (5.75" W x 2" L x 2" H)



**Safepak® 900 (formerly Safepak® XL)**

ACC-9140

Holds up to 5 standard cryovial boxes (5" W x 5" L x 2" H). Holds up to 15 Cryoport cryovial boxes (5.75" W x 2" L x 2" H)



**Cryostrap®**

ACC-9141

Holds up to 9 cassettes (250 ml). Holds up to 12 Cryoport cryovial boxes (5.75" W x 2" L x 2" H). Holds up to 5 standard cryovial boxes (5" W x 5" L x 2" H)

**Tamper Evidence & Security** - Please insert quantity to order



**Serialized Seal**

ACC-9140

Serialized blue tamper-evident seal that helps identify if shipment has been compromised.



**Cryovial Boxes**

ACC-9043

Holds up to 30, 1.2 or 2.0 ml cryovials



**Standard Canister**

ACC-9149

Holds up to 14 specialty canes



**Safepak® Soft System 1800**

SET-1000

Includes Safepak® 1800 Bag, High Volume Sleeve, and SoftRack™



**Safepak® System 1800**

SET-1001

Includes Safepak® 1800 Bag and High Volume Sleeve



**Cable Zip Tie**

ACC-10631

Serialized blue tamper-evident cable zip ties that helps identify if shipment has been compromised.

### Additional Information:

By completing and submitting this form to Cryoport Systems, LLC, you acknowledge and agree that you have read and agree to the **Terms and Conditions** set forth at [www.cryoport.com/legal](http://www.cryoport.com/legal).